



H

The Arc Wi 2G adapter is expanding the functionality of Arc sensors by providing wireless communication for local monitoring in parallel to robust 4-20 mA signal, and simple sensor connection to the PCS with additional internal galvanic isolator for an enhanced signal quality.

This mobile device empowers the operator to monitor measurement values, calibrate Arc sensors and configure various parameters with a unified user interface for pH, DO, Conductivity and ORP. The Arc View Mobile device is based on the Samsung Galaxy Tab Active tablet and comes pre-configured with the ArcAir application, app blocker application, power supply cable, instruction manual and Hamilton quick guide.

	Description
1	Arc View Mobile Basic
3	Arc View Mobile Advanced

The Deutschmann Unigate[®] CL Module converts the Modbus protocol of the Arc and VisiFerm DO sensors into the Profibus DP protocol. All Arc parameters DO, pH, ORP and conductivity are supported. The conversion script is pre-installed.

	Ref
Profibus Converter	243555
Profibus Programmer's Manual	624719

The Arc Modbus OPC converter converts the Modbus protocol of Arc and VisiFerm DO sensors into OPC UA protocol. All Arc parameters are