

# Ordering information continued.

For drive units with E-coupling

—

Sales No.	Product	ø / length mm	Volume ml	max. Tip speed m/s
<b>X «X-DESIGN» - FOR HARD &amp; BRITTLE PILLS AND GEL CAPS</b>				
11030255	PT-DA 20 / 2XEC-E116	20 / 116	min. 10*	21
11030277	PT-DA 20 / 2XEC-E157	20 / 157	min. 10*	21
*When using a vessel ø 22 mm / volume is always dependent on the diameter of the vessel				
<b>Z «Z-DESIGN» - FOR PILLS AND HARD / BRITTLE SAMPLES</b>				
11030067	PT-DA 12 / ZMEC-E157*	12 / 157	2 - 250	12
* Z-Design aggregate with knives				
<b>B «BIOTRONA® -DESIGN» - HIGH-TURBULENCE MIXING</b>				
11030034	PT-DA 07 / BEC-E107	07 / 107	10 - 100	7
11030286	PT-DA 12 / BEC-E157	12 / 157	50 - 500	12
11030288	PT-DA 20 / BEC-E192	20 / 192	20 - 2000	21
11030291	PT-DA 25 / BEC-E192	25 / 192	250 - 4000	26
<b>TiN «TITANIUM NITRATED» - UPGRADE</b>				
11035050	TiN rotor for the selected dispersing aggregate			ø 12 - 60 mm
11035150	TiN stator for the selected dispersing aggregate			ø 12 - 60 mm
<b>G «GAS TIGHT - WITH SINGLE MECHANICAL SEAL» FOR PRESSURE / VACUUM</b>				
11030273	PT-DA 20 / 2G-E210	20 / 210	20 - 2000	22