




Funnel Adaptor	
 <p>190100-74</p>	Adaptor for Sartorius 16A07 plastic funnel
 <p>190100-75</p>	Adaptor for Pall 4870 (100ml), 4871 (250ml) filter funnel
 <p>190100-76</p>	Hose barb adaptor (for Ø4-8mm tubing)

Maintenance

1. This product is not autoclavable. Please clean the case by pure water or 75% alcohol.
2. It is strongly recommended to replace O-ring of funnel adaptor once every year at least to keep good aspiration.
3. Silicone tubes are consumables; please replace at least every half year to ensure good drainage.