

## designed for scientists



## **MMT 100.1**

/// Data Sheet

The Multi-use Milling Tube MMT 100.1 with a maximum volume of 100 ml can be used and cleaned in a dishwasher multiple times. The package includes a rich set of spare parts so that wearing parts can be replaced, if necessary. Therefore it is possible to carry out a number of experiments with the grinding chamber depending on the nature of the sample.

Scope of delivery: 1 x MMT 100, 25 x sealing, 5 x beater, 5 x coupling

Additionally, there is an abrasion set A-MMT.100 available.











## designed for scientists

## **Technical Data**

Weight [kg] 0.2816







