



designed for scientists



## S 10 N - 10 G - ST Dispersing tool

/// Data Sheet

Suitable for dispersing instrument T 10 basic and T 10 standard

ST tools are recommended for use with tissue and other fibrous materials\*

Working range: 1 - 100 ml

Stator diameter: 10 mm

Immersion depth: 22 - 100 mm

Materials in contact with medium: PTFE, AISI 316L

Suitable for solvents

Max. temperature: 180 °C

Sterilization methods: all methods



designed for scientists

### Technical Data

Volume range (H <sub>2</sub> O) [l]	0.001 - 0.1
Stator diameter [mm]	10
Rotor diameter [mm]	7.5
Gap between rotor and stator [mm]	0.25
allowable Speed max. [rpm]	30000
Circumferential speed max. [m/s]	11.8
Immersion depth [mm]	23 - 100
Shaft length [mm]	115
Material in contact with medium	PTFE, AISI 316L
pH range	2 - 13
Suitable for solvents	yes
Suitable for abrasive substances	yes
Working temperature max. [°C]	180
Sterilization methods	all methods
Ultimate fineness, suspensions [µm]	5 - 25
Ultimate fineness, emulsions [µm]	1 - 10
Weight [kg]	0.075
Permissible ambient temperature [°C]	5 - 40