

Swing-out Rotor, GRS-G-r750-4

4 loadings, 90° | for 1580, 1580R

Bucket, GLB-r750-r750

3,600 rpm / 2,898 xg




Bucket, GLB-500-r750

3,600 rpm / 3,159 xg



Microplate Holder, GLP-mw-r750

3,600 rpm / 2,485 xg



Tube Racks and Adaptors for GLB-r750-r750



GAM-5-30(r750)
30 x 3/5 ml



GAM-15-20(r750)
20 x 10/15 ml



GAM-c15-16(r750)
16 x 15 ml (C)



GAM-50-8(r750)
8 x 50 ml



GAM-c50-5(r750)
5 x 50 ml (C)



GAS-250(r750)
250 ml



GAS-c250(r750)
250 ml (C)



GAS-c500(r750)
500 ml (C)

Swing-out Rotor, GRS-G-750-4

4 loadings, 90° | for 1580, 1580R

Bucket, GLB-750-750

4,000 rpm / 3,667 xg



Bucket, GLB-b750-750 with O-ring sealed Cap

4,000 rpm / 3,667 xg




Bucket, GLB-500-750

4,000 rpm / 3,900 xg



Microplate Holder, GLP-mw-750

4,000 rpm / 3,068 xg



Tube Racks and Adaptors for GLB-750-750 / GLB-b750-750



GAM-5-24(r750)
24 x 3/5 ml



GAM-10-21(r750)
21 x 10 ml



GAM-15-19(r750)
19 x 15 ml



GAM-c15-14(r750)
14 x 15 ml (C)



GAM-50-7(r750)
7 x 50 ml



GAM-c50-5(r750)
5 x 50 ml (C)




GAS-250(r750)
250 ml









GAS-c250(r750)
250 ml (C)






GAS-c500(r750)
500 ml (C)



Ordering Information - Adaptors and Multi-Tube Racks for Swing-Out Rotors

Cat. No. & Description		Capacity	Hole Dimension (ϕ xL, mm)	Applicable Rotor & Bucket
Adaptors and Multi-tube Racks	with O-ring sealed Cap			
GAS-c500(r750) Adaptor for 500 ml conical	-	1 x 500 ml conical	98 x 94	GLB-r750-r750 
GAS-c250(r750) Adaptor for 250 ml conical	-	1 x 250 ml conical	61 x 93	
GAS-250(r750) Adaptor for 250 ml	-	1 x 250 ml	62.5 x 87	
GAM-c50-5(r750) Multi-tube rack for 50 ml conical	-	5 x 50 ml conical	29.5 x 85	
GAM-50-8(r750) Multi-tube rack for 50 ml	-	8 x 50 ml	29.5 x 85	
GAM-c15-16(r750) Multi-tube rack for 15 ml conical	-	16 x 15 ml conical	17.5 x 85	
GAM-15-20(r750) Multi-tube rack for 15 ml	-	20 x 10/15 ml	17.2 x 85	
GAM-5-30(r750) Multi-tube rack for 3/5 ml	-	30 x 3/5 ml	13 x 56.5 (3 ml) 13 x 85 (5 ml)	GLB-750-750  GLB-b750-750 
GAS-c500(750) Adaptor for 500 ml conical	-	1 x 500 ml conical	95.6 x 58	
GAS-c250(750) Adaptor for 250 ml conical	-	1 x 250 ml conical	60.5 x 125	
GAS-250(750) Adaptor for 250 ml	OK	1 x 250 ml	62.3 x 87	
GAM-c50-5(750) Multi-tube rack for 50 ml conical	OK	5 x 50 ml conical	29.2 x 88	
GAM-50-7(750) Multi-tube rack for 50 ml	OK	7 x 50 ml	29.2 x 87	
GAM-c15-14(750) Multi-tube rack for 15 ml conical	OK	14 x 15 ml conical	17 x 85	
GAM-15-19(750) Multi-tube rack for 15 ml	OK	19 x 15 ml	17.2 x 85	
GAM-10-21(750) Multi-tube rack for 10 ml	OK	21 x 10 ml	16 x 85	
GAM-5-24(750) Multi-tube rack for 3/5 ml	OK	24 x 3/5 ml	13.2 x 58 (3 ml) 13.2 x 85 (5 ml)	
GAS-250(r250) Adaptor for 250 ml	-	1 x 250 ml	62.5 x 87	
GAM-85-2(r250) Multi-tube rack for 85 ml & 15 ml	-	2 x 85 & 15 ml	17 x 86.5 (15 ml) 38.5 x 86.5 (85 ml)	
GAM-c50-3(r250) Multi-tube rack for 50 ml conical	-	3 x 50 ml conical	30.5 x 90	
GAM-50-4(r250) Multi-tube rack for 50 ml	-	4 x 50 ml	30.2 x 90	
GAM-s50-4(r250) Multi-tube rack for 50 ml skirted	-	4 x 50 ml skirted	31.4 x 65	
GAM-c15-9(r250) Multi-tube rack for 15 ml conical	-	9 x 15 ml conical	17.5 x 90	
GAM-15-12(r250) Multi-tube rack for 10/15 ml	-	12 x 10/15 ml	17.5 x 90	
GAM-5-12(r250) Multi-tube rack for 3/5 ml	-	12 x 3/5 ml	13.5 x 58 (3 ml) 13.5 x 80 (5 ml)	
GAM-m2.0-20(r250) Multi-tube rack for 1.5/2.0 ml	-	20 x 1.5/2.0 ml	11.3 x 39	GLB-60-250 GLB-d20/10-250 
GAS-60s(60) Adaptor for 60 ml syringe	-	1 x 60 ml syringe	Custom Size	
GAS-20s(20) Adaptor for 20 ml syringe	-	1 x 20 ml syringe		
GAS-10s(10) Adaptor for 10 ml syringe	-	1 x 10 ml syringe		
GAS-c50(50) Adaptor for 50 ml conical	OK	1 x 50 ml conical	29.5 x 15.6	GLB-b50-100 
GAS-u15(50) Adaptor for 15 ml round & conical	OK	1 x 15 ml round & conical	17.3 x 87	
GAS-100(b100) Adaptor for 100 ml	OK	1 x 100 ml	44.2 x 93	GLB-b100-100 
GAS-85(b100) Adaptor for 85 ml	OK	1 x 85 ml	38.5 x 96.4	
GAS-c50(b100) Adaptor for 50 ml conical	OK	1 x 50 ml conical	30 x 100	
GAS-50(b100) Adaptor for 50 ml	OK	1 x 50 ml	29.5 x 95.9	
GAS-c15(b100) Adaptor for 15 ml conical	OK	1 x 15 ml conical	17.2 x 106.5	
GAM-c15-3(100) Multi-tube rack for 15 ml conical	-	3 x 15 ml conical	17.2 x 106.5	
GAM-15-3(100) Multi-tube rack for 9~15 ml	-	3 x 9 ~ 15 ml	17.2 x 105	
GAM-10-5(b100) Multi-tube rack for 4~10 ml	OK	5 x 4 ~ 10 ml	16 x 60	
GAM-7-5(b100) Multi-tube rack for 2~7 ml	OK	5 x 2.6~7 ml	13.5 x 60	
GAM-m2.0-6(b100) Multi-tube rack for microtube	OK	6 x 1.5/2.0 ml	11 x 39	

Catalog Numbering Rule of Rotors and Accessories



Type Symbol of Tube or Bucket

- No symbol : round type
- r : rectangular type
- c : conical type
- d : dual hole type
- b : O-ring sealed bucket
- u : universal for both of round & conical tube
- mw : for multiwell plate