
































 50 ml	 50 ml conical	 85 ml	 250 ml	 500 ml	 1,000 ml		
6		6		6		6	
GRF-SR-50-6 ① GRF-L-50-6 ②	GRF-G-50-8 ① GRF-L-50-8 ②	GRF-SR-c50-6 ① GRF-G-c50-6 ② GRF-L-c50-6 ③	GRF-G-85-6 ① GRF-L-85-6 ②	GRF-G-250-6 ① GRF-L-250-6 ②	GRF-L-500-6	GRF-L-1000-4	GRF-L-1000-6
							
29.5 120	29.5 120	29.5 125	38.5 122 ① / 125 ②	62 (250 ml) / 171(5 ml) 130 (250 ml) / 120 (15 ml)	70 170	98 175	98 175
-	GAS-15(50) 15 ml GAS-c15(60) 15 ml conical	GAS-c15(50) 15 ml conical	GAS-15(85) 15 ml GAS-c15(85) 15 ml conical GAS-50(85) 50 ml GAS-c50(85) 50 ml conical	-	GAS-250(500) 250 ml	GAS-250(1000) 250 ml GAS-500(1000) 250 ml	GAS-250(1000) 250 ml GAS-500(1000) 250 ml
6,000 / 3,824 ①	-	6,000 / 3,960 ①	-	-	-	-	-
-	12,000 / 16,582 ①	10,000 / 11,101 ②	10,000 / 12,298 ①	-	-	-	-
-	15,000 / 25,910 ①	13,000 / 18,592 ②	12,000 / 17,709 ①	8,000 / 9,588 ①	8,000 / 11,305	-	-
17,000 / 32,310 ②	15,000 / 25,910 ②	15,000 / 24,753 ②	15,000 / 25,029 ②	10,000 / 14,981 ②	10,000 / 14,981 ②	8,000 / 11,305	-
22,000 / 54,111 ③	20,000 / 46,062 ③	17,000 / 31,793 ③	20,000 / 44,496 ③	13,000 / 25,318 ③	10,000 / 17,664	8,000 / 12,071	7,000 / 10,825
-	-	-	-	10,000 / 14,981 ③	8,000 / 11,305	-	-

- **A** : Adaptor available for smaller tubes
- (C) : conical tube
- The fixed angle rotors are supplied with rotor lid (except PCR strip rotors).
- The open style angle rotors are supplied with tube sleeves.
- The swing-out rotor may be ordered together with compatible buckets and adaptors.
- The colors of the rotors, buckets and adaptors may be different and changeable without notification.

 50 ml	 50 ml conical	 100 ml	 100 ml glass	 250 ml	 250 ml conical	 400 ml blood bag	 500 ml	 500 ml conical	 750 ml
1 x 4		1 x 4 or 1 x 6		1 x 4		1 x 4		1 x 4	
GLB-50-50	GLB-b50-100	GLB-c50-50	GLB-b100-100 (with GAS-100(100)) GLB-100(L)-250 GLB-100(C)-250	GLB-250-250 ① GLB-c250-250 ② (with GAS-250(250))	GLB-750-750 ① ② GLB-b750-750 ③ ④ (with GAS-250(750))	GLB-c250-250	GLB-750-750	GLB-1400-1400 (with GAM-400-2(1400))	GLB-500-750
GRS-S-50-4 ① GRS-SR-50-4 ②	GRS-S-100-4 ① GRS-SR-100-4 ② GRS-G-100-6	GRS-S-50-4 ① GRS-SR-100-4 ② GRS-G-100-6	GRS-S-100-4 ① GRS-SR-100-4 ② GRS-G-100-6	GRS-G-250-4 ① GRS-G-r250-4 ② GRS-L-250-4 ③	GRS-G-750-4 ① GRS-L-750-4 ② GRS-L-750-6 ③	GRS-G-r750-4 ① GRS-L-r750-4 ②	GRS-G-750-4 ① GRS-L-750-4 ② GRS-L-750-6 ③	GRS-L-1400-4W GRS-L-750-4 ① GRS-L-750-4 ② GRS-L-750-6 ③	GRS-G-750-4 ① GRS-L-750-4 ② GRS-L-750-6 ③
									
-	-	125	120/130 (w/o cap) 207 (L) 190(S, C)	153 ① / 130 ② ③	130 / 142(w/o cap)	135	153	150	165
GAS-u15(50) 15 ml & 15 ml (c) GAS-c50(50) 50 ml conical	-	-	GAM-m2-0(250) 2.0 ml GAM-7-5(100) 7 ml GAM-10-5(100) 10 ml GAM-15-3(100) 15 ml (c) GAM-c15-3(100) 15 ml (c) GAS-50(100) 50 ml (c) GAS-c50(100) 50 ml (c) GAS-85(100) 85 ml GAS-100(100) 100 ml	GAM-m2-0(250) 2.0 ml GAM-5-12(250) 5 ml GAM-15-12(250) 15 ml (c) GAM-c15-9(250) 15 ml (c) GAM-50-4(250) 50 ml (c) GAM-c50-3(250) 50 ml (c) GAM-85-2(250) 85 ml GAS-250(250) 250 ml	GAM-5-24(750) 5 ml GAM-10-21(750) 10 ml GAM-15-19(750) 15 ml (c) GAM-c15-14(750) 15 ml (c) GAM-50-7(750) 50 ml (c) GAM-c50-5(750) 50 ml (c) GAS-250(750) 250 ml GAS-c250(750) 250 ml (c) GAS-500(750) 500 ml (c)	GAM-5-30(750) 5 ml GAM-15-20(750) 15 ml (c) GAM-c15-16(750) 15 ml (c) GAM-50-8(750) 50 ml (c) GAM-c50-5(750) 50 ml (c) GAS-250(750) 250 ml GAS-c250(750) 250 ml (c) GAS-500(750) 500 ml (c)	GAM-400-2(1400) 400 ml GAM-200-3(1400) 200 ml double bag GAM-5-30(750) 5 ml GAM-15-20(750) 15 ml (c) GAM-c15-16(750) 15 ml (c) GAM-50-8(750) 50 ml (c) GAM-c50-5(750) 50 ml (c) GAS-250(750) 250 ml GAS-c250(750) 250 ml (c) GAS-500(750) 500 ml (c)	GAM-5-24(750) 5 ml GAM-10-21(750) 10 ml GAM-15-19(750) 15 ml (c) GAM-c15-14(750) 15 ml (c) GAM-50-7(750) 50 ml (c) GAM-c50-5(750) 50 ml (c) GAS-250(750) 250 ml (c) GAS-c250(750) 250 ml (c) GAS-500(750) 500 ml (c)	GAM-5-24(750) 5 ml GAM-10-21(750) 10 ml GAM-15-19(750) 15 ml (c) GAM-c15-14(750) 15 ml (c) GAM-50-7(750) 50 ml (c) GAM-c50-5(750) 50 ml (c) GAS-250(750) 250 ml (c) GAS-c250(750) 250 ml (c) GAS-500(750) 500 ml (c)
4,000 / 2,700 ① 4,800 / 3,888 ②	4,000 / 2,700 ① 4,200 / 2,976 ②	4,000 / 2,700 ① 4,800 / 3,888 ②	4,000 / 2,774 ① 4,200 / 3,059 ②	-	-	-	-	-	-
-	4,000 / 3,032 ③	-	4,000 / 3,105 ③	5,000 / 5,394 ① 4,000 / 3,134 ②	-	5,000 / 5,394	-	-	-
-	-	-	-	4,000 / 3,667 ①	3,600 / 2,898 ①	4,000 / 3,667 ①	3,600 / 2,898 ①	4,000 / 3,900 ① 4,000 / 4,248 ②	3,600 / 3,159 ①
-	-	-	-	-	-	-	-	4,000 / 3,667 ① 4,000 / 3,982 ②	3,600 / 2,898 ①
-	-	-	4,000 / 4,358 (L) 4,000 / 4,155 (S) 4,000 / 4,191 (C)	4,000 / 3,667 ② 4,000 / 3,982 ③	3,600 / 2,898 ②	-	4,000 / 3,982	4,000 / 3,982 ② 4,000 / 3,982 ③	3,600 / 2,898 ② 4,000 / 3,667 ②
-	-	-	-	4,000 / 3,667 ③	3,600 / 2,898 ③	-	4,000 / 3,982	4,000 / 3,982 ③	4,000 / 3,667 ③ 4,000 / 3,982 ③