



# CHAPTER 15 PCR Consumables


Our tube design and thickness offers the perfect seal, since the seal between the tube opening and the tube cap is critical in minimizing evaporation or sample leakage. Loose caps allow evaporation; overly tight caps cause spillage upon opening. All PCR tubes are non-sterile, certified RNase and DNase free.

Available in a wide variety of styles from individual tubes, to strips and plates.  
Most of PCR plates have both clear and white type.  
Numeric markings on our strip of tubes.  
Alphanumeric markings on plates.



### Non-skirt, qPCR plates, unskirted, pp

Part NO.	Option	Color	Pack
PCRP-20NS	0.2ml 96-well	clear	10/box, 10boxes/case
PCRP-10NS	0.1ml 96-well	clear	10/box, 10boxes/case
PCRP-20NS-W	0.2ml 96-well	white	10/box, 10boxes/case
PCRP-10NS-W	0.1ml 96-well	white	10/box, 10boxes/case




### Half skirt, fit ABI, qPCR plates, pp

Part NO.	Option	Color	Pack
PCRP-20HS	0.2ml 96-well	clear	10/box, 10boxes/case
PCRP-20HS-W	0.2ml 96-well	white	10/box, 10boxes/case




### Half skirt, fit Biorad/Roche, qPCR plates, pp

Part NO.	Option	Color	Pack
PCRP-10HSR	0.1ml 96-well	clear	10/box, 10boxes/case
PCRP-10HSR-W	0.1ml 96-well	white	10/box, 10boxes/case




### Sub semi-skirt or High-skirt, fit ABI, qPCR plates, pp

Part NO.	Option	Color	Pack
PCRP-20IS	0.2ml 96-well	clear	10/box, 10boxes/case
PCRP-10IS	0.1ml 96-well	clear	10/box, 10boxes/case




### Full-skirt, qPCR plates, pp

Part NO.	Option	Color	Pack
PCRP-10FS	0.1ml 96-well	white base, clear wells	10/box, 10boxes/case
PCRP-10FS-W	0.1ml 96-well	white base, white wells	10/box, 10boxes/case
PCRP-20FS	0.2ml 96-well	white base, clear wells	10/box, 10boxes/case
PCRP-20FS-W	0.2ml 96-well	white base, white wells	10/box, 10boxes/case



40ul 384-Well PCR plates, pp

Part NO.	Color	Pack
PCRP-40FS-384C	white base, clear wells	10/box, 10boxes/case
PCRP-40FS-384C-W	white base, white wells	10/box, 10boxes/case



PCR sealing film (works for all kinds of plates), individual pack

Part NO.	Option	Color	Pack
PCRF-CP	Advanced sealing film for qPCR	clear	100/unit



0.2ml 8-Strip PCR Tubes, with Clear Cap (qPCR), pp

Part NO.	Color	Cap	Pack
PCRS20F	clear tubes	flat cap	125/unit, 10units/case
PCRS-20F-W	white tubes	flat cap	125/unit, 10units/case
PCRS-20D	clear tubes	round cap	125/unit, 10units/case



0.1ml 8-Strip PCR Tubes, with Clear Cap (qPCR), pp

Part NO.	Color	Cap	Pack
PCRS-10F	clear tubes	flat cap	125/unit, 10units/case
PCRS-10F-W	white tubes	flat cap	125/unit, 10units/case




0.2ml PCR Tubes with Attached Flat Caps, pp

Part NO.	Color	Pack
PCRT-20F	clear tubes	1,000/unit, 10units/case

# CHAPTER 16

# Compatibility List (PCR Consumable)







PART NO.		Roche					Abbott	Agilent/Stratagena				Analytik Jena/Biometra						Bioer								Xian Tianlong Science				TaKaRa Bio Inc				LongGene																					
		Cobas z480	LightCycler®480 System 384-well	LightCycler®480 System	LightCycler®96 System	LightCycler®Nano	m2000rt	Mx3000P™	Mx3005P™	Mx4000™	Agilent AriaMX	SureCycler 8800 thermal cycler	Toptical	qTower	Tone96/96G	Alpha SC 0.2ml*96 0.2ml*48 0.2ml*48*2	FlexCycler 0.2ml*96 □ 0.2ml*48	TRIO48 0.2ml*48	TRIO Combi 0.2ml*48	SpeedCycler 36	Line-Gene I/Mini 16	Line-Gene K48	Line-Gene 3310/3320	Line-Gene 9600	Line-Gene 9600plus	GenePro B	Life ECO Thermal Cycler	Life Express TC-96/T/H(a) □ TC-48/H(t)	LifeProTC-96/G/H(b)A	Life Touch TC-96/G/H(b)B	XPcycler TC-XP-	Gentier 96E	DTC-4	DTC-4G	TL988 Real-time PCR	Thermal Cycler Dice TP600	Thermal Cycler Dice TP650	Thermal Cycler Dice TP850	Thermal Cycler Dice TP870	A100	A200	Gradient PCR	L Series 96	Multi-Purpose PCR	MyGene™ PCR	Q2000C									
PCR 8 StripTube	PCRS-10F									√					√	√										√	√												√	√	√	√	√	√	√	√									
	PCRS-10F-W	√	√	√	√	√				√		√				√	√																																√						
	PCRS-20F						√	√	√	√	√		√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√						
	PCRS-20F-W						√	√	√	√	√		√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√						
PCR Tube	PCRT-20F																																																						
PCR Plate	PCRP-20NS									√						√	√										√	√																						√	√	√	√	√	√
	PCRP-20NS-W	√	√	√	√	√				√		√				√	√																																				√		
	PCRP-10NS						√	√	√	√	√		√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√		
	PCRP-10NS-W						√	√	√	√	√		√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√	√		
	PCRP-10HSR	√r	√r	√r	√r	√r																																																	
	PCRP-10HSR-W	√r	√r	√r	√r	√r																																																	
	PCRP-20HS							√																																															
	PCRP-20HS-W							√																																															
	PCRP-20IS																																																						
	PCRP-10IS																																																						
PCRP-10FS																																																							

Note: √r r=recommend