

compact ROCKER™

- An examples to use
- Western blotting
 - Gel staining
 - Immunoprecipitation
 - DNA / RNA Extraction

■ CR95 / Economical and compact rocker

- Speed setting of 12, 15, 21, 27, 39 or 55rpm
- Suitable for continuous operation
- Can be used in incubators up to 80°C

■ CR100 / Compact rocker in timed and continuous operation

- Variable speed control from 8 and 70rpm
- Continuous or timed operation

■ CR300 / Large capacity platforms and smooth rocking motion

- Variable speed control from 8 and 70rpm
- Continuous or timed operation
- Large platform for variety of flasks and tubes

► Common facts

- Sturdy and durable construction for use at variable temperature
- Non-slip mat as standard, optional dimpled pad prevents tubes from rolling around
- All moving parts of the rocking platform are combined with ball joints for smooth operation
- Optional second platform available to double the working area



► CR300



► CR100



► CR95




Specification

Model	CR95	CR100	CR300
Tilt Angle	±7°		±9°
Speed Range	selected 12 ~ 55rpm (6steps)		8 ~ 70 rpm
Motor	DC Warm Gear Motor		
Timer (speed)	-		180min
Temp. Resist	Lower than 80°C(Dry air)		Lower than 65°C(Dry air)
Loadable Weight	5Kg	7Kg	10Kg
Weight	7Kg	5Kg	7Kg
Dimensions (mm)	W310xD200xH160	W310xD200xH170	W330xD300xH203
Platform (mm)	W290xD200	W290xD200	W300xD300
Wattage	23W		35W
Power	AC 100-240V (*220-230V : CE cert)		Free Voltage (AC100-240V)
Frequency	50/60Hz		50/60Hz

Accessories

Model	Model Number	Description
CR95	DP-95	Dimpled Pad (W290xD200 mm)
	TY95	Second platform (with posts : 4ea)
CR100	DP-100	Dimpled Pad (W290xD200 mm)
	TY100	Second platform (with posts : 4ea)
CR300	DP-300	Dimpled Pad (W300xD300 mm)
	TY300	Second platform (with posts : 4ea)
	UP400	Upper Platform (W400xD400 mm)

ACCESSORIES (CR95)




DP-95



TY95



ACCESSORIES (CR100)




DP-100



TY100



ACCESSORIES (CR300)



DP-300



TY300



UP400