

### Additional Rotors (Swing Out):



|                                  |
|----------------------------------|
| <b>Item: Z306-100</b>            |
| 4 x 145ml Swing Out Rotor        |
| Max Speed: 4,500 rpm / 3,350 x g |
| Accel: 15sec. / Decel: 15 sec.   |

| Inserts:       | Tube            | qty/<br>insert | qty./<br>rotor |
|----------------|-----------------|----------------|----------------|
| Z306-100-A145  | 145ml           | 1              | 4              |
| Z306-100-A100: | 100ml           | 1              | 4              |
| Z306-100-A85:  | 85ml            | 1              | 4              |
| Z306-100-A50:  | 50ml<br>conical | 1              | 4              |
| Z306-100-A15:  | 15ml<br>conical | 2              | 8              |
| Z306-100-A10:  | 10ml            | 7              | 28             |
| Z306-100-A57:  | 5/7ml           | 7              | 28             |
| Z306-100-A20:  | 1.5/2.0ml       | 10             | 40             |

\*Inserts sold in 2/pk.  
Rotor accepts 4 inserts (2 x 2/pk.)



|  |
|--|
| <b>Item: Z326-100H</b>   |
| 4 x 100ml Swing Out Rotor <b>with</b><br>hermetically sealed buckets |
| Max Speed: 5,000 rpm / 3,885 x g                                     |
| Accel: 13 sec. / Decel: 13 sec.                                      |

| Inserts         | Tube            | qty/<br>insert | qty./<br>rotor |
|-----------------|-----------------|----------------|----------------|
| Z326-100H-A100: | 100ml           | 1              | 4              |
| Z326-100H-A85:  | 85ml            | 1              | 4              |
| Z326-100H-A50:  | 50ml<br>conical | 1              | 4              |
| Z326-100H-A15:  | 15ml<br>conical | 2              | 8              |
| Z326-100H-A10:  | 10ml            | 4              | 16             |
| Z326-100H-A57:  | 5/7ml           | 4              | 16             |
| Z326-100H-A20:  | 1.5/2.0ml       | 5              | 20             |

\*Inserts for Z326-100H Rotor, sold in 2/pk.  
Rotor accepts 4 inserts (2 x 2/pk.)

|  |
|--|
| <b>Item: Z326-100H-AC</b>                |
| Same as above,<br><b>without</b> buckets |



| Buckets        | Tube  | qty/<br>Bucket | qty./<br>rotor |
|----------------|---|----------------|----------------|
| Z326-100H-B50: | 50ml  | 2              | 8              |
| Z326-100H-B15: | 15ml  | 3              | 12             |
| Z326-100H-BUK: | Hermetically Sealed Buckets (for use<br>with Z326-100H-XXX inserts above) |                |                |

\*Buckets for Z326-100H-AC rotor, sold in 2/pk.  
Rotor accepts 4 buckets (2 x 2/pk.)

### Also Available from Benchmark:

## Z36 HK

### Super Speed Refrigerated Centrifuge

- Convenient benchtop design
- Produces g-forces up to 65,395 x g
- 10 versatile rotor options

