

Accessories Z326 and Z326-K

Rotor/Bucket

Order No.:
Z306-100



Swing out rotor 4 place
Angle: 0-90°

Adapter

Acceleration in s: 15/100
Deceleration in s: 15/150
max. Speed: 4500 rpm
RCF-value: 3350 xg
max. Radius: 14.8 cm
Temp. at max. speed: -1°C

Volume
Tubes ø in mm

145 ml Nalgene® 50	100 ml RB Glass 44	85 ml RB/OR 38	50 ml RB Glass 29	50 ml Falcon® 29
--------------------------	--------------------------	----------------------	-------------------------	------------------------



Length min./max. in mm
Tubes per rotor
Radius correction (cm)
Order-No.: Z306-100-

100 4 1.0 A145	90/1 4 1.0 A100	90/108 4 0.8 A85	90/118 4 0.7 A50R	116 4 0.3 A50
-------------------------	--------------------------	---------------------------	----------------------------	------------------------

Volume

Tubes ø in mm

50 ml skirted 29	30 ml Universal 25	15 ml Falcon® 17	6-10 ml Vakut.® 16	1.8-7 ml Vakut.® 13	1.5 ml Reakt. 11
------------------------	--------------------------	------------------------	--------------------------	---------------------------	------------------------



Length min./max. (mm)
Tubes per rotor
Radius Correction (cm)
Order No. Tube rack

116 4 0.6 A50SS	95/105 4 1.1 A30	120 8 0.0 A15	90/108 8 0.4 A10	75/100 28 0.7 A57	42 40 5.5 A20
--------------------------	---------------------------	------------------------	---------------------------	----------------------------	------------------------

Order No.:
Z326-100H

Swing out rotor 4 Places
w/ sealed buckets
Angle: 0-90°

Z326-100H-BUK
Bucket for Tube Rack w/ lid

Z326-100H-LID
Lid for bucket



Acceleration in s: 15/100
Deceleration in s: 15/150
max. Speed: 5000 rpm
RCF-value: 4080 xg
max. Radius: 14.6 cm
Temp. at max. speed: 2°C

Volume
Tubes ø in mm

100 ml RB Glass 40	85 ml RB/OR 38	50 ml RB/OR 29	50 ml Falcon® 29
--------------------------	----------------------	----------------------	------------------------



Length min./max. in mm
Tubes per rotor
Radius correction (cm)
Order-No.: Z326-100H-

90/108 4 0.9 A100	90/108 4 0.6 A85	90/133 4 0.6 A50R	116 4 0.4 A50
----------------------------	---------------------------	----------------------------	------------------------

Volume

Tubes ø in mm

50 ml skirted 29	30 ml Universal 25	15 ml Falcon® 17	15 ml Falcon® 17	15 ml RB Glass 17	6-10 ml Vakut.® 16	10-15ml Vakut.® 16	1.8-7 ml Vakut.® 13	1.5 ml Reakt. 11	50ml Falcon® 29	15ml Falcon® 17
------------------------	--------------------------	------------------------	------------------------	-------------------------	--------------------------	--------------------------	---------------------------	------------------------	-----------------------	-----------------------



Length min./max. (mm)
Tubes per rotor
Radius Correction (cm)
Order No. Tube rack

116 4 0.7 A50SS	95/105 4 1.5 A30	120 4 0.5 A15S	120 8 0.5 A15*	90/108 16 0.6 A15R	75/100 16 0.6 A10	75/125 4 0.6 A15L*	75/100 16 0.6 A57	42 20 0.6 A20	116 8 0.0 B50	120 12 0.0 B15
--------------------------	---------------------------	-------------------------	-------------------------	-----------------------------	----------------------------	-----------------------------	----------------------------	------------------------	------------------------	-------------------------

* only without bucket lid

Order No.:
Z326-02MP



Swing out rotor 2 x 3 Microplates
Angle: 0-90°

Acceleration in s: 20/160
Deceleration in s: 15/360
max. Speed: 4500 rpm
RCF-value: 2716 xg
max. Radius: 12.0 cm
max. Loading height: 63 mm
max. Loading weight: 310 g
Temp. at max. speed: 0°C

Order No.:
BT1500-A1520
BT1500-A05
Z326-MPCA

Tube rack
20 x 1.5 ml (2.0 ml) Tubes
20 x 0.5 ml Tubes
3 x Chamber slide or 35/60 mm petri dish

Order No.:
Z326-3015



Angle rotor 30 x 15 ml
Angle: 35°

Acceleration in s: 15/160
Deceleration in s: 10/380
max. Speed: 4500 rpm
RCF-value: 2830/2467 xg
max. Radius: 12.5/10.9 cm
max. Ø: 17 mm
Temp. at max. speed: -6°C

Order No.:
without
C3100-ADP

Tube rack	Ø (mm)	Length (min./max.):
1 x 15 ml-Glass	17.0	95/115
1 x 5/7 ml (pack of 8)	13	75/100