

Accessories Z326 and Z326-K

Rotor/Bucket

Order No.:
Z326-6405



Angle rotor 64 x 0.5 ml
Angle: 40°

Adapter

Acceleration in s: 17/150
Deceleration in s: 12/170
Max. speed: 13,500 rpm
Max. RCF: 16,298/17,725 xg
Max. radius: 8.0/8.7 cm
Max. Ø: 8 mm
Temp. at max. speed: 6°C

Order No.:
Z326-3202



Angle rotor 4 x 8 PCR Strips
Angle: 45°

Acceleration in s: 16/130
Deceleration in s: 12/130
Max. speed: 15,000 rpm
Max. RCF: 15,343 xg
Max. radius: 6.1 cm
Max. Ø: 6 mm
Temp. at max. speed: 5°C

Order No.:
Z326-1205



Angle rotor 12 x 5 ml Reaction Tubes
Angle: 45°C

Acceleration in s: 25/205
Deceleration in s: 16/223
Acceleration in s: 27/220
Deceleration in s: 17/373
Max. speed: 14,000 rpm
Max. speed: 15,000 rpm
Max. RCF: 18,624 xg
Max. RCF: 21,379 xg
Max. radius: 8.5 cm
Max. ø: 17 mm
Temp. at max. speed: 6°C

Order No.:

Z216-05-C10 1.0 ml Cryo tubes (ø 12,6 x 41 mm)
set of adapters (6 pcs.)
Z216-05-C18 1.8 ml Cryo tubes (ø 12,6 x 48 mm)
set of adapters (6 pcs.)

The intuitive control panel makes it easy for anyone to operate **Hermle Centrifuges** straight away.

