

# Accessories Z327 and Z327-K

## Rotor/Bucket

Order No.:  
**Z327-200**

Swing out rotor 4 place w/ buckets  
Angle: 0-90°  
(w/o buckets: Z327-200-AC)

Z327-200-BUK  
Bucket for Adapter

Z327-200-LID  
Lid for Bucket | hermetically sealed



## Adapter

Acceleration in s: 21/102  
Deceleration in s: 12/387  
max. speed: 5,000 rpm  
RCF-Value: 4,136 xg  
max. radius: 14.8 cm

Volume

Tubes Ø in mm

200 ml  
bottle  
56

200 ml  
conical  
56

100 ml  
RB Glass  
40

85 ml  
RB/OR  
38



Length min./max. (mm)

Tubes per Rotor

Radius Correction (cm)

Order No. Adapter: Z327-200-\_\_\_

112  
4  
0.0  
without

120  
4  
0.0  
A200

95/115  
4  
0.2  
A100

95/108  
4  
0.2  
A85

Volume

Tubes Ø in mm

50 ml  
RB/OR  
29

50 ml  
conical  
29

50 ml  
skirted  
29

25/30 ml  
Universal  
25

25/30 ml  
RB/OR  
25

15 ml  
conical  
17

15 ml  
conical  
17

15 ml  
RB  
17

5/7 ml  
RB  
13

5 ml  
conical  
17

1.5/2.0 ml  
micro  
11

6-10 ml  
Vakut.®  
16

6-10 ml  
Vakut.®  
16

1.8-7 ml  
Vakut.®  
13

50 ml  
conical  
29



Length min./max. (mm)

Tubes per rotor

Radius Correction (cm)

Order No. Adapter: Z327-200-\_\_\_

95/103  
4  
0.2  
A50R

116  
4  
0.1  
A50

116  
4  
0.2  
A50SS

95/105  
4  
0.2  
A25

90/110  
4  
0.2  
A30

120  
16  
0.1  
A15B\*

120  
12  
0.0  
A15

60/107  
28  
0.2  
A15R

65/100  
32  
0.2  
A57

60/70  
20  
0.1  
A5F

40  
32  
0.1  
A20

130  
4  
0.2  
A15L

100  
28  
0.2  
A10

75/100  
32  
0.2  
A57

116  
8  
0.0  
B50

Order No.:  
**Z306-100**



Swing out rotor 4 place  
Angle: 0-90°

Acceleration in s: 12/93  
Deceleration in s: 12/226  
max. speed: 4,500 rpm  
RCF-Value: 3,350 xg  
max. radius: 14.8 cm

Volume

Tubes Ø in mm

145 ml  
Nalgene®  
50

100 ml  
RB Glass  
44

85 ml  
RB/OR  
38

50 ml  
RB Glass  
29

50 ml  
conical  
29



Length min./max. (mm)

Tubes per rotor

Radius Correction (cm)

Order No. Adapter: Z306-100-\_\_\_

100  
4  
1.0  
A145

90/108  
4  
1.0  
A100

90/108  
4  
0.8  
A85

90/118  
4  
0.7  
A50R

116  
4  
0.3  
A50

Volume

Tubes Ø in mm

50 ml  
Rim  
29

30 ml  
Universal  
25

15 ml  
conical  
17

15 ml  
RB  
17

6-10 ml  
Vakut.®  
16

1.8-7 ml  
Vakut.®  
13

1.5/2.0 ml  
micro  
11



Length min./max. (mm)

Tubes per rotor

Radius Correction (cm)

Order No. Adapter

116  
4  
0.6  
A50SS

95/105  
4  
1.1  
A30

120  
8  
0.0  
A15

90/108  
28  
0.4  
A10

90/108  
28  
0.4  
A10

75/100  
28  
0.7  
A57

40  
40  
5.5  
A20

Order No.:  
**Z326-100H**

Swing out rotor 4 Places  
w/ sealed buckets  
Angle: 0-90°

Z326-100H-BUK  
Bucket for Tube Rack w/ lid

Z326-100H-LID  
Lid for bucket



Acceleration in s: 15/100  
Deceleration in s: 15/150  
max. Speed: 5000 rpm  
RCF-value: 4080 xg  
max. Radius: 14.6 cm  
Temp. at max. speed: 2°C

Volume

Tubes ø in mm

100 ml  
RB Glass  
40

85 ml  
RB/OR  
38

50 ml  
RB/OR  
29

50 ml  
Falcon®  
29



Length min./max. in mm

Tubes per rotor

Radius correction (cm)

Order-No.: Z326-100H-\_\_\_

90/108  
4  
0.9  
A100

90/108  
4  
0.6  
A85

90/133  
4  
0.6  
A50R

116  
4  
0.4  
A50

Volume

Tubes Ø in mm

50 ml  
skirted  
29

30 ml  
Universal  
25

15 ml  
Falcon®  
17

15 ml  
Falcon®  
17

15 ml  
RB Glass  
17

6-10 ml  
Vakut.®  
16

10-15 ml  
Vakut.®  
16

1.8-7 ml  
Vakut.®  
13

1.5 ml  
Reakt.  
11

50 ml  
Falcon®  
29

15 ml  
Falcon®  
17



Length min./max. (mm)

Tubes per rotor

Radius Correction (cm)

Order No. Tube rack

116  
4  
0.7  
A50SS

95/105  
4  
1.5  
A30

120  
4  
0.5  
A15S

120  
8  
0.5  
A15\*

90/108  
16  
0.6  
A15R

75/100  
16  
0.6  
A10

75/125  
4  
0.6  
A15L\*

75/100  
16  
0.6  
A57

42  
20  
0.6  
A20

116  
8  
0.0  
B50

120  
12  
0.0  
B15

\* only without bucket lid