



Accessories Z327 and Z327-K

Rotor/ Bucket		Adapter					
Order No.: Z327-200 Swing out rotor 4 place w/ buckets Angle: 0-90° (w/o buckets: Z327-200-AC) Z327-200-BUK Bucket for Adapter Z327-200-LID Lid for Bucket hermetically sealed		Acceleration in s: 21/102 Deceleration in s: 12/387 max. speed: 5,000 rpm RCF-Value: 4,136 xg max. radius: 14.8 cm					
		Volume Tubes Ø in mm		200 ml bottle 56	200 ml conical 56	100 ml RB Glass 40	85 ml RB/OR 38
		Length min./max. (mm) Tubes per Rotor Radius Correction (cm) Order No. Adapter: Z327-200-___		112 4 0.0 without	120 4 0.0 A200	95/115 4 0.2 A100	95/108 4 0.2 A85

Volume	50 ml RB/OR 29	50 ml conical 29	50 ml skirted 29	25/30 ml Universal 25	25/30 ml RB/OR 25	15 ml conical 17	15 ml conical 17	15 ml RB 17	5/7 ml RB 13	5 ml conical 17	1.5/2.0 ml micro 11	6-10 ml Vakut.® 16	6-10 ml Vakut.® 16	1.8-7 ml Vakut.® 13	50 ml conical 29
Tubes Ø in mm															
Length min./max. (mm)	95/103	116	116	95/105	90/110	120	120	60/107	65/100	60/70	40	130	100	75/100	116
Tubes per rotor	4	4	4	4	4	16	12	28	32	20	32	4	28	32	8
Radius Correction (cm)	0.2	0.1	0.2	0.2	0.2	0.1	0.0	0.2	0.2	0.1	0.1	0.2	0.2	0.2	0.0
Order No. Adapter: Z327-200-___	A50R	A50	A50SS	A25	A30	A15B*	A15	A15R	A57	A5F	A20	A15L	A10	A57	B50

Order No.:		Adapter						
Z306-100 Swing out rotor 4 place Angle: 0-90°		Acceleration in s: 12/93 Deceleration in s: 12/226 max. speed: 4,500 rpm RCF-Value: 3,350 xg max. radius: 14.8 cm						
		Volume Tubes Ø in mm		145 ml Nalgene® 50	100 ml RB Glass 44	85 ml RB/OR 38	50 ml RB Glass 29	50 ml conical 29
		Length min./max. (mm) Tubes per rotor Radius Correction (cm) Order No. Adapter: Z306-100-___		100 4 1.0 A145	90/108 4 1.0 A100	90/108 4 0.8 A85	90/118 4 0.7 A50R	116 4 0.3 A50

Volume	50 ml Rim 29	30 ml Universal 25	15 ml conical 17	15 ml RB 17	6-10 ml Vakut.® 16	1.8-7 ml Vakut.® 13	1.5/2.0 ml micro 11
Tubes Ø in mm							
Length min./max. (mm)	116	95/105	120	90/108	90/108	75/100	40
Tubes per rotor	4	4	8	28	28	28	40
Radius Correction (cm)	0.6	1.1	0.0	0.4	0.4	0.7	5.5
Order No. Adapter	A50SS	A30	A15	A10	A10	A57	A20

Order No.:		Adapter					
Z326-100H Swing out rotor 4 Places w/ sealed buckets Angle: 0-90° Z326-100H-BUK Bucket for Tube Rack w/ lid Z326-100H-LID Lid for bucket		Acceleration in s: 15/100 Deceleration in s: 15/150 max. Speed: 5000 rpm RCF-value: 4080 xg max. Radius: 14.6 cm Temp. at max. speed: 2°C					
		Volume Tubes ø in mm		100 ml RB Glass 40	85 ml RB/OR 38	50 ml RB/OR 29	50 ml Falcon® 29
		Length min./max. in mm Tubes per rotor Radius correction (cm) Order-No.: Z326-100H-___		90/108 4 0.9 A100	90/108 4 0.6 A85	90/133 4 0.6 A50R	116 4 0.4 A50

Volume	50 ml skirted 29	30 ml Universal 25	15 ml Falcon® 17	15 ml Falcon® 17	15 ml RB Glass 17	6-10 ml Vakut.® 16	10-15 ml Vakut.® 16	1.8-7 ml Vakut.® 13	1.5 ml Reakt. 11	50 ml Falcon® 29	15 ml Falcon® 17
Tubes Ø in mm											
Length min./max. (mm)	116	95/105	120	120	90/108	75/100	75/125	75/100	42	116	120
Tubes per rotor	4	4	4	8	16	16	4	16	20	8	12
Radius Correction (cm)	0.7	1.5	0.5	0.5	0.6	0.6	0.6	0.6	0.6	0.0	0.0
Order No. Tube rack	A50SS	A30	A15S	A15*	A15R	A10	A15L*	A57	A20	B50	B15