

Accessories Z446 and Z446-K

Rotor/Bucket

Order No.:
Z446-750-AC*
Swing out rotor 4 place
Angle: 0-90°



*Z446-750
Swing out rotor with
round buckets included

Z446-750-LID*
PC-Screw cap
for round bucket



Adapter

Order No.:
Z446-750-BUK
Round bucket



Acceleration in s: 47/483
Deceleration in s: 49/1,287
Acceleration in s: 55/551
Deceleration in s: 54/1,501
max. speed: 4,000/4,500 rpm
max. RCF: 3,434/4,346 xg
max. radius: 19.2 cm
Temp. at max. speed: 4°C

Z446-MP-LID*
PC-Lid
for rectangular bucket



Z446-MP-RA
Removing aid
for rectangular bucket




Order No.:
Z446-MP-BUK
Rectangular bucket
for Microplates



Acceleration in s: 47/483
Deceleration in s: 49/1,287
Acceleration in s: 55/551
Deceleration in s: 54/1,501
max. speed: 4,000/4,500 rpm
max. RCF: 3,130/3,961 xg
max. radius: 17.5 cm
max. loading height: 80 mm
max. loading weight: 640 g
Temp. at max. speed: 4°C

Adapters for rectangular bucket
on request

Order No.:
Z446-750-MP
Microplate carrier



Acceleration in s: 47/482
Deceleration in s: 49/1,287
Acceleration ins s: 55/551
Deceleration in s: 54/1,501
max. speed: 4,000/4,500 rpm
mac. RCF: 2,754/3,486 xg
max. radius: 15.4 cm
max. loading height: 72 mm
max. loading weight: 310 g
Temp. at max. speed: -3°C

*hermetically sealable

Volume	Corning®				Corning®	Falcon®	conical	conical	50 ml	30 ml	conical	conical	15 ml
	750 ml	500 ml	500 ml	250/290 ml	250 ml	175/225 ml	50 ml	50 ml			15 ml	15 ml	
ø in mm	99	99	97	62	60	61	29	29	29	26	17	17	17
Round bucket													
Length min./max. (mm)	152	114	148	122/137	159	118/137	116	116	102/123	102/123	120	120	97/120
Tubes per rotor	4	4	4	4	4	4	20	28	28	32	56	68 ¹	72
Radius correction (cm)	0	0	0.5	0.2	0.4	0.2	0.2	0.2	0.2	0.2	0.1	0.3	0.5
Order No. Adapter Z446-750	without	without	-A500	-A250	-A250F	-A225	-A50-E	-A50	-A50R	-A30	-A15-E	-A15	-A15R
¹ 2-level													
	Vacutainer®												
Volumen	1.8-5 ml	4-7 ml	5 ml	6-10 ml									
ø in mm	13	13	16	16									
Round bucket													
Length min./max. (mm)	75	100	75	100									
Tubes per rotor	84	84	72	72									
Radius correction (cm)	0.3	0.3	0.3	0.3									
Order No. Adapter Z446-750	-A57	-A57	-A10	-A10									

Order No.:
Z326-02MP



Swing out rotor 2 x 3 Microplates
Angle: 0-90°

Acceleration in s: 27/251
Deceleration in s: 23/602
max. speed: 4,500 rpm
max. RCF: 2,716 xg
max. radius: 12 cm
max. loading height: 63 mm
max. loading weight: 310 g
Temp. at max. speed: -8°C

Order No.:
BT1500-A1520
BT1500-A05
Z326-MPCA

Tube rack
20 x 1.5 ml (2.0 ml) Tubes
20 x 0.5 ml Tubes
3 x Chamber slide or 35/60 mm petri dish

Order No.:
Z326-3020H



Angle rotor 30 x 1.5/2.0 ml
Angle: 45°

Acceleration in s: 24/218
Deceleration in s: 30/526
max. speed: 15,000 rpm
max. RCF: 23,643 xg
max. radius: 9.4 cm
max. ø: 11 mm
Temp. at max. speed: 1°C

Order No.:	Tube rack	ø (mm)	Length (min./max.)
without	1 x 2,0 ml	11	40
without	1 x 1,5 ml	11	40
C1008-A5-6	1 x 0.5 ml	8	45
C1008-A2-6	1 x 0.4/0.2 ml PCR	6	45