



designed for scientists



## MMT 40.1

/// Data Sheet

Multi-use Milling Tube MMT 40.1

The Multi-use Milling Tube MMT40.1 with a maximum volume of 40 ml can be used and cleaned in a dishwasher multiple times. The package includes a rich set of spare parts so that wearing parts can be replaced if necessary. Therefore it is possible to carry out a number of experiments with the grinding chamber depending on the nature of the sample.

Delivery: 1 x MMT 40, 25 x sealing, 5 x beater, 5 x coupling

[www.ika.com](http://www.ika.com)

Subject to technical changes



IKAworlwide



IKAworlwide /// #lookattheblue



@IKAworlwide



designed for scientists

Additional there is an Abrasion set A-MMT.100 available.





designed for scientists

**Technical Data**

---

Weight [kg]

0.262

---

