

# Ordering Information

## Foxy R2 and Foxy R1

Foxy R1 fraction collector with diverter valve: Includes 1/16-inch OD (0.031 ID) tubing for collection up to 25 mL/min and (1) rack for 13 mm tubes. Other racks may be substituted; specify racks when ordering .....68-3870-012

Foxy R2 fraction collector with diverter valve: Includes 1/8-inch OD (0.062 ID) tubing for collection up to 125 mL/min and (2) racks for 13 mm tubes. Other racks may be substituted; specify racks when ordering .....68-2130-009

Foxy R2 fraction collector with high-flow diverter valve: Includes compression fittings for 1/4-inch OD tubing. Maximum flow rate 1 L/min. Includes (2) racks for 25 mm tubes. Other racks may be substituted; specify rack when ordering .....68-2130-007

### Extra Racks for Foxy R1 & Foxy R2 (Order 2 racks for Foxy R2.)

13 mm diameter tubes (144).....	60-2137-150
16 mm diameter tubes (100).....	60-2137-151
18 mm diameter tubes (72).....	60-2137-152
25 mm diameter tubes (36).....	60-2137-153
1.5 ml microcentrifuge tubes (60).....	60-2137-154
29 mm centrifuge tubes (36).....	60-2137-157
28 mm diameter scintillation vials (36).....	60-2137-156
96-well microplates (2).....	60-2137-159

Set of two racks for (18) 480 mL glass bottles. Each rack holds (9) bottles. Fits Foxy R2 only; requires qty. 1 set .....60-5247-004

96-well microplates (6)  
For Foxy R2 only, requires qty. 1 .....60-2137-139

Preparative funnel rack with 26 ETFE\* or ECTFE\*\* funnels and 50 ft. of 3/8-inch ID vinyl, FEP-lined tubing. For Foxy R2 only, requires qty. 1. Organic compatible .....60-2137-162

Preparative funnel rack with 36 ETFE\* funnels and 50 ft. of 1/4-inch ID vinyl tubing. (Requires 1 rack for Foxy R1; 2 racks for Foxy R2). Non-organic.....60-2137-163

### Interface Cables for Foxy R1 & Foxy R2


HP1100/1200 system for remote start signal and peak separation.....60-5394-474

Non-Teledyne ISCO pump connect cable .....69-2134-170

External control of Foxy R1 or R2 by contact closure or TTL I/O.....69-2134-171

### Accessories for Foxy R1 & Foxy R2

Tubing and fitting package for 25–125 mL/min  
Glass bottles, 480 mL, case of 40 .....60-2137-062



For complete information about Teledyne ISCO chromatography products, visit;

[www.teledyneisco.com](http://www.teledyneisco.com)

\* Equivalent to Tefzel, a registered trademark of E.I. DuPont de Nemours & Co.

\*\* Equivalent to Halar, a registered trademark of Allied Chemical Corp.