


Accessories for CLG3-70 / PSR3-70

LC GreenBox

- Real-time monitoring of equipment connected to LC Connected possible through the dedicated mobile application for Lab Companion's.
(Supported from Android 6.0 Marshmallow version and above. iOS devices to be supported in the future)
- Each lab member can access connected devices for which they have access permission through their personal mobile app.
- One LC GreenBox can communicate with up to four Lab Companion's devices via an RS-232 port and can be expanded via multiple LC GreenBoxes.
- Real-time alarms for events concerning connected equipment.
(temperature deviation, long-term door opening, power failure, and etc.)



LC GreenBox



Real-time monitoring through mobile devices



Connection up to 4 devices



Power failure notification

Cat. No.	Description	Dimension (W x D x H, mm / inch)
AAAQ1011	LC GreenBox	156 x 94 x 34 / 6.1 x 3.7 x 1.3

Battery

- Activate/deactivate battery via settings.
- When the setting of battery is ON, a warning alarm occurs if the power is turned off during operation.
- The recorder also records power failure events.

Cat. No.	Description
AAAE15301	Battery