

# Accessories Z287-A

## Rotor/Bucket

Order No.:  
Z287-0450H



Swing out rotor 4 x 50 ml  
Angle: 0-75°

Order No.:  
Z287-0430H



Swing out rotor 4 x 30 ml Universal  
Angle: 0-90°

Order No.:  
Z287-0810H



Swing out rotor 8 x 10 ml  
for blood collection tubes  
Angle: 0-90°

Order No.:  
Z287-02MP



Swing out rotor 2 x Microtitre plates  
Angle: 0-90°

Order No.:  
Z287-0650



Angle rotor 6 x 50 ml  
Angle: 40°

Order No.:  
Z287-1215



Angle rotor 12 x 15 ml  
Angle: 32°

## Adapter

Acceleration in s: 26/242  
Deceleration in s: 33/252  
max. speed: 5,000 rpm  
RCF-value: 3,633 xg  
max. radius: 13.0 cm  
max. Ø: 29 mm

| Order No.:      | Tube rack         |
|-----------------|-------------------|
| without adapter | 1 x 50 ml conical |
| Z326-50-B15     | 1 x 15 ml conical |
| Z326-15-A7*     | 1 x 7 ml          |
| Z326-15-A5*     | 1 x 5 ml          |

| Ø [mm] | Length (min./max.): |
|--------|---------------------|
| 29     | 116                 |
| 17     | 120                 |
| 13.5   | 95/115              |
| 13.5   | 65/90               |

\* only possible with adapter Z326-50-B15

Acceleration in s: 26/241  
Deceleration in s: 31/250  
max. speed: 5,000 rpm  
RCF-value: 3,437 xg  
max. radius: 12.3 cm  
max. Ø: 25 mm

| Order No.:      | Tube rack         |
|-----------------|-------------------|
| without adapter | 1 x 30 ml skirted |
| without adapter | 1 x 25 ml skirted |

| Ø [mm] | Length (min./max.): |
|--------|---------------------|
| 25     | 95                  |
| 25     | 90                  |

Acceleration in s: 25/240  
Deceleration in s: 32/253  
max. speed: 5,000 rpm  
RCF-value: 3,717 xg  
max. radius: 13.3 cm  
max. Ø: 16 mm

| Order No.:   | Tube rack                |
|--------------|--------------------------|
| Without      | 1 x 6-10 ml Vacutainer®  |
| Z326-0810-A7 | 1 x 1.8-7 ml Vacutainer® |

| Ø [mm] | Length (min./max.): |
|--------|---------------------|
| 16     | 100                 |
| 13     | 75/100              |

Acceleration in s: 27/182  
Deceleration in s: 29/193  
max. speed: 3,800 rpm  
RCF-value: 1,453 xg  
max. radius: 9.0 cm  
max. loading height: 20 mm  
max. loading weight: 100 g

Acceleration in s: 34/293  
Deceleration in s: 38/317  
max. speed: 6,000 rpm  
RCF-value: 4,427 xg  
max. radius: 11.0 cm  
max. Ø: 29 mm

| Order No.:      | Tube rack         |
|-----------------|-------------------|
| without adapter | 1 x 50 ml conical |
| without adapter | 1 x 50 ml         |
| Z206-50-A15     | 1 x 15 ml conical |
| Z206-50-A57     | 1 x 7 ml          |
| Z206-50-A57C*   | 1 x 5 ml          |

| Ø [mm] | Length (min./max.): |
|--------|---------------------|
| 29     | 116                 |
| 29     | 100/120             |
| 17     | 120                 |
| 13     | 95/115              |
| 13     | 65/90               |

\* With Z206-50-A57

Acceleration in s.: 34/295  
Deceleration in s: 38/317  
max. speed: 6,000 rpm  
RCF-value: 4,427 xg  
max. radius: 11.0 cm  
max. Ø: 17 mm

| Order No.:      | Tube rack         |
|-----------------|-------------------|
| without adapter | 1 x 15 ml conical |
| without adapter | 1 x 15 ml glass   |
| Z326-15-A7      | 1 x 7 ml          |
| Z326-15-A5      | 1 x 5 ml          |

| Ø [mm] | Length (min./max.): |
|--------|---------------------|
| 17     | 120                 |
| 17     | 105/120             |
| 13     | 95/115              |
| 13     | 65/90               |